

Active protection with high security

# SERAFIT® PROTECT

Surgical sutures

Novel sutures  
with antibac-  
terial coating



150  
YEARS  
1866 - 2016



SERAG  
WIESSNER



## SERAFIT® PROTECT

Antibacterial coating

Good knot stability

Outstanding suppleness


Minimum sawing effect

SERAFIT® PROTECT is an absorbable braided suture that is outstandingly supple, passes smoothly through the tissues, but has a minimal saw effect. It can be tied reliably, even in areas containing blood and secretions, and ensures good knot stability.

### Antibacterial coating

SERAFIT® PROTECT may help to prevent infection. The possibility of bacterial colonisation of the thread is considerably reduced. The antibacterial coating consists of chlorhexidine, a broad spectrum antibiotic, which has been used in the cosmetic industry and medicine for many years.

### Results of agar diffusion tests






*Staphylococcus aureus*



*Methicillin-resistant Staphylococcus epidermidis*



Material	 POLYGLYCOLIC ACID
Symbol	 violet, multifilament (braided), coated
Strength	USP 3/0 bis 1 EP 2 bis 4
Absorption profile	50% tensile strength after 15-20 days no residual tensile strength after 60-90 days
Available combinations	without needle: single thread / multipacks with needle: with DS, FRX, HR, HRT, HRX, HS
Indications / areas of use	ligatures / dermatology / gastroenterology / gynaecology / urology




In vitro absorption profile  
SERAFIT® PROTECT violet  
EP 2, USP 3/0


**150**  
YEARS  
1866 - 2016




**SERAG**  
**WIESSNER**


SERAG-WIESSNER GmbH & Co. KG  
Zum Kugelfang 8 - 12  
95119 Naila/Germany


 + 49 9282 937-0


 + 49 9282 937-9369

*Export Department:*

 + 49 9282 937-230

 + 49 9282 937-9785

 [info@serag-wiessner.de](mailto:info@serag-wiessner.de)

 [www.serag-wiessner.de](http://www.serag-wiessner.de)